



香港中文大學
The Chinese University of Hong Kong



香港中文大學醫學院
Faculty of Medicine
The Chinese University of Hong Kong



Oroeasy Manual for Trainer



Supported by Knowledge Transfer Fund at CUHK
香港中文大學知識轉移基金支持項目

Download Oroeasy for Mobile Device



Oroeasy

The Department of Otorhinolaryngology, Head and Neck Surgery of the Chinese University of Hong Kong (CUHK) is launched an Oroeasy mobile app in February 2021. This mobile app is funded by the Social Innovation and

apps.apple.com



Oroeasy - Apps on Google Play

The Department of Otorhinolaryngology, Head and Neck Surgery of the Chinese University of Hong Kong (CUHK) is launched an Oroeasy mobile app in February 2021. This mobile app is funded by the Social Innovation and

play.google.com



v1.0.3 (build 9) on iTunes App Store

<https://apps.apple.com/us/app/oroeasy/id1560136173>

v1.0.3 (ver. code 9) on Google Play

<https://play.google.com/store/apps/details?id=hk.edu.cuhk.ent.oroeasy>

Outline

Login as a Trainer

Assign the Training
to Trainee

Create a Trainee

Create a Training
Template



香港中文大學
The Chinese University of Hong Kong



香港中文大學醫學院
Faculty of Medicine
The Chinese University of Hong Kong




Login as a Trainer

When the app is launched, select “Trainer”



Trainee



Trainer

English

Change Language here

QrOeasy



香港中文大學
The Chinese University of Hong Kong



香港中文大學醫學院
Faculty of Medicine
The Chinese University of Hong Kong



Login as trainer by your registered email and password

Trainer

Login

Please provide your information.



Email

xxxxxx@xxx.xxx



Password

XXXXXXXXXXXX

OR

Create an account

Oroeasy app contains comprehensive video database for oral motor exercises. Registered trainer could create individualised treatment regime for their patients. Speech therapists and their appointee are eligible to register the account. For more information, please visit the Oroeasy official webpage.

Sign In

Create an account

English

Oroeasy



香港中文大學
The Chinese University of Hong Kong



香港中文大學醫學院
Faculty of Medicine
The Chinese University of Hong Kong



40th Anniversary
1978-2018

If you have not registered for the trainer account, click “create an account” that direct you to the official webpage for registration.

Trainer

Login

Please provide your information.



Email



Password

Sign In

OR

Create an account

Oroeasy app contains comprehensive video database for oral motor exercises. Registered trainer could create individualised treatment regime for their patients. Speech therapists and their appointee are eligible to register the account. For more information, please visit the Oroeasy official webpage.



oroeasy.com

Create an account



香港中文大學
The Chinese University of Hong Kong



香港中文大學醫學院
Faculty of Medicine
The Chinese University of Hong Kong



40th Anniversary
1978-2018

If you forgot the password, click “forgot password” to retrieve it.

Trainer

Login

Please provide your information.



Email

XXXXXX@XXX.XXX



Password

....

OR

Create an account

Oreasy app contains comprehensive video database for oral motor exercises. Registered trainer could create individualised treatment regime for their patients. Speech therapists and their appointee are eligible to register the account. For more information, please visit the Oreasy official webpage.

Sign In

Create an account

Forgot password?

English

Oreasy



香港中文大學
The Chinese University of Hong Kong



香港中文大學醫學院
Faculty of Medicine
The Chinese University of Hong Kong



40th Anniversary
1978-2018

By entering your registered email, you will receive an email to prompt you to setup a new password.

Retrieve your password

Please enter a valid email address.

Email

Email XXXXXX@XXX.XXX

Submit



香港中文大學
The Chinese University of Hong Kong




香港中文大學醫學院
Faculty of Medicine
The Chinese University of Hong Kong



Create a Training Template

Click “  ” to access the menu panel,
then select “Training Collection”



Welcome XXX

Exercise Lips

Training Collection

- Contain the training template

Trainee List

- Contain the trainees that are under your care

Language

- Switch among 3 language
i.e. English/Traditional Chinese/Simplified Chinese

Logout

- return to the main page

Pause Time No. of Repetition



Click “ + Training Set” create a training

Create Training Set

(en-us)New - 2021-04-19 10:22:30

Cancel

Confirm

+ Training Set



(en-us)New - 2021-04-19 10:22:30

Duration 0 min No. of Exercise 0

Category

Description
 1. (en-us)Description goes here.

- Delete
- Duplicate
- Edit

Click "Edit" to edit the training.

Training Set

Enter Training name and training description.

Details

TrainingSet name:

Lips Strengthening

To strengthen the muscles of the lips

訓練集名稱:

強化嘴唇

強化嘴唇肌肉

训练集名称:

強化咀唇

強化嘴唇肌肉

Duration
0 min

Break Time Between Exercises

5 s ×



Update Training Set

1. Please make sure that you have input all 3 languages – English, traditional Chinese and simplified Chinese.
2. If all of your clients use only one language, you could put the boxes in the same language.

Set the break time between exercises.
Click “Update Training Set” to save the input.

Details

TrainingSet name:

Lips Strengthening

To strengthen the muscles of the lips

訓練集名稱:

強化嘴唇

強化嘴唇肌肉

训练集名称:

強化咀唇

強化嘴唇肌肉

 Duration
0 min

Break Time Between Exercises

5 s ▾ ✕



Update Training Set



Click “+ Exercise” to insert new exercise to training

Details

TrainingSet name:

Lips Strengthening

To strengthen the muscles of the lips

訓練集名稱:

強化嘴唇

強化嘴唇肌肉

训练集名称:

強化咀唇

強化嘴唇肌肉

 Duration
0 min

Break Time Between Exercises

5 s  



 Update Training Set

No exercise in this training set.

 Exercise

Change Exercise Category

The screenshot shows the 'Add Exercise' screen in an application. At the top, there is a blue header bar with a hamburger menu icon on the left, the text 'Add Exercise' in the center, and a dropdown menu on the right currently set to 'Cheek'. A red circle highlights the 'Cheek' dropdown, and a purple arrow points to it from a dark grey menu that is open on the right side of the screen. This menu lists several categories: Cheek, Face, Jaw, Lips, Soft Palate, Speech, and Tongue. Below the header, there are two exercise cards. The first card is for 'Cheek compress (ID: 101)' and the second is for 'Cheek stretch (ID: 102)'. Each card displays details such as Action Time (10 s), Pause Time (5 s), No. of Repetition (0), and Equipment (none). The 'Cheek' category is also listed for both exercises. Below the details, there are sections for 'Procedures' with numbered steps. At the bottom left of the screen, there is a green button labeled 'Exercise' with a list icon.

Add Exercise Cheek

Cheek compress (ID: 101)

Action Time	10 s	Pause Time	5 s	No. of Repetition	0
Category	Cheek	Equipment	none		

Procedures

1. Close your mouth
2. Suck in cheeks and hold

Cheek stretch (ID: 102)

Action Time	10 s	Pause Time	5 s	No. of Repetition	0
Category	Cheek	Equipment	none		

Procedures

1. Take a breath
2. Tightly close your lips
3. Puff the air full in your mouth and hold

Exercise

Cheek
Face
Jaw
Lips
Soft Palate
Speech
Tongue



Click “+ Exercise” to add exercise to training



Add Exercise

Lips



Lip close (ID: 401)

 Action Time
10 s

 Pause Time
5 s

 No. of Repetition
0

 Category
Lips

 Equipment
none

Procedures

 1. Press your lips tightly together and hold

 Exercise

Lip close with resistance (ID: 402)

 Action Time
10 s

 Pause Time
5 s

 No. of Repetition
0

 Category
Lips

 Equipment
Tongue depressor



Select model, viewing angle and no. of repetition

Lip close (ID: 401)

Action Time

Pause Time

No. of Repetition

Add Exercise

Viewing Angle and Model

 Viewing Angle and Model ▾ ✕

No. of Repetition

 Select no. of Repetition ▾ ✕

Cancel

Confirm



After “confirm”, the exercise is added to training.
Note: it may take time to refresh the training from server

TrainingSet name:

Lips Strengthening

To strengthen the muscles of the lips

訓練集名稱:

強化嘴唇

強化嘴唇肌肉

训练集名称:

強化咀唇

強化嘴唇肌肉

Duration

1 min

Break Time Between Exercises

5 s ▾ ✕



Update Training Set

Lip close (ID: 401)

Action Time
10 s

Pause Time
5 s

No. of Repetition
4

Viewing Angle
45°

Model
Male Adult

Category
Lips

Equipment
none

Exercise

Procedures



Click “+ Exercise” to insert more exercise to training

Details

TrainingSet name:

Lips Strengthening

To strengthen the muscles of the lips

訓練集名稱:

強化嘴唇

強化嘴唇肌肉

训练集名称:

強化咀唇

強化嘴唇肌肉

 Duration
0 min

Break Time Between Exercises

5 s  



 Update Training Set

No exercise in this training set.

 + Exercise

Edit exercises in the training template

☰ Training Set Info



✎ Lip close (ID: 401)

🏃 Action Time
10 s

⏸ Pause Time
5 s

👁 Viewing Angle
45°

👤 Model
Male Adult

📌 Category
Lips

🔧 Equipment
none

Procedures

❓ 1. Press your lips tightly together and hold

Edit - to edit the 'Viewing Angle and Model'
Up - to arrange the order of the exercise (move up)
Down - to arrange the order of the exercise (move down)
Delete - to delete the exercise

📺 Preview

✎ Edit

↑ Up

↓ Down

🗑 Delete

✎ Lip close with resistance (ID: 402)

🏃 Action Time
10 s

⏸ Pause Time
5 s

🔄 No. of Repetition
4

👁 Viewing Angle
45°

👤 Model
Male Adult

📌 Category
Lips

🔧 Equipment
Tonque depressor

Note: it may take time to refresh the training from server





香港中文大學
The Chinese University of Hong Kong




香港中文大學醫學院
Faculty of Medicine
The Chinese University of Hong Kong



Create a Trainee

Click “  ” to access the menu panel,
then select “Trainee List”



Welcome XXX

Exercise Lips

Training Collection

Trainee List

Language

About Oroeasy

App Policy

Logout

Pause Time No. of Repetition

Training Collection

- Contain the training template

Trainee List

- Contain the trainees that are under your care

Language

- Switch among 3 language
i.e. English/Traditional Chinese/Simplified Chinese

Logout

- return to the main page



The page lists the trainees who registered under your account.
Click “+ Trainee” to create a trainee

☰ Trainee List Active ▾

📷 Reference Code AAA

Registration Status
In Progress

👤 Detail ☰ Training Set

+ Trainee

We Care About the Privacy of your Patient!

Create Trainee

Password

 Please enter a 4- to 8-digit Password. If leave empt...

Name

 (Optional) Trainee Name Goes Here

Reference Code

 (Mandatory) Trainee Reference Code Goes Here

Gender

 Select Gender ▼ ✕

Year of Birth

 Select Year of Birth ▼ ✕

Diagnosis - Primary

 (Optional) Trainee Primary Diagnosis Goes Here

Diagnosis - Secondary

 (Optional) Trainee Secondary Diagnosis Goes Here

+ Trainee

- Only Reference code is mandatory in this form which helps you identify your client, so the code should be unique in your account.
- The other info is optional.

When filled in the form, click “+ Trainee” to create the account.

Create Trainee

Password



12345678

Name



CHAN TT

Reference Code



CUHK 001

Gender



MALE ▾ ×

Year of Birth



1931 ▾ ×

Diagnosis - Primary



(Optional) Trainee Primary Diagnosis Goes Here

Diagnosis - Secondary



(Optional) Trainee Secondary Diagnosis Goes Here

 Trainee



Click “Detail” for login information

◀ Search 2:30 PM Mon 26 Apr 23% 🔋

☰ Trainee List Active ▾

🖼️ Reference Code CUHK 001

Registration Status
Completed

👤 Detail 📄 Training Set

🖼️ Reference Code AAA

Registration Status
In Progress

👤 Detail 📄 Training Set

Trainee requires trainee ID and password to login to their account.
Note: the *leading zero* in trainee ID is not required

Details

 Trainee ID
00000005

Password
 12345678

 Name
CHAN TT

 Gender
MALE

Diagnosis - Primary
 -

Diagnosis - Secondary
 -

 Update Trainee Info

Training ID: 5
Password: 12345678

It is suggested to provide the Trainee Login Form with these info to the client.

The form could be downloaded from the official webpage Oroeasy.com. The form includes:

1. A quick guide
2. The QR code of Oroeasy app store
3. Place for filling in trainee ID and password



香港中文大學
The Chinese University of Hong Kong




香港中文大學醫學院
Faculty of Medicine
The Chinese University of Hong Kong



Assign the Training to Trainee

Click “  ” to access the menu panel,
then select “Training Collection”



Welcome XXX

Exercise Lips

Training Collection

- Contain the training template

Trainee List

- Contain the trainees that are under your care

Language

- Switch among 3 language
i.e. English/Traditional Chinese/Simplified Chinese

Logout

- return to the main page

Pause Time No. of Repetition



Choose the training template to be assigned to the patient, then click “Edit”



Training Collection

Lips Strengthening

 Duration
3 min

 No. of Exercise
3

 Category
Lips

Description

 1. To strengthen the muscles of the lips


 Delete

 Duplicate

 Edit

Click “  ” for training assignment

☰ Training Set Info



 Details

TrainingSet name:

Lips Strengthening

To strengthen the muscles of the lips

訓練集名稱:

強化嘴唇

強化嘴唇肌肉

训练集名称:

強化咀唇

強化嘴唇肌肉

 Duration
3 min

Break Time Between Exercises

 5 s 

 Update Training Set

 Lip close (ID: 401)

 Action Time
10 s

 Pause Time
5 s

 No. of Repetition
4

Viewing Angle

Model

Select Trainee ID, Training Start Date, End Date and No. of Trial

TrainingSet name:

訓練集名稱:

训练集名称:

Lips Stre

To streng

Durati

3 min

Break

5 s

Update

Assign Training Set to Trainee

(id: 00000003) AAA ▾ ×



Trial Target

3 ▾ ×



Start Date

2021-04-22 ×



End Date

2021-04-30 ×



Cancel

Confirm



Click “confirm” to assign the training to the trainee

Lips Strengthening	強化嘴唇	強化咀唇
To strengthen the muscles of the lips	強化嘴唇肌肉	強化嘴唇肌肉

Duration
3 min

Break Time Between Exercises
5 s

Assign Training Set to Trainee
Done

Update Training Set

Lip close (ID: 401)

Action Time 10 s	Pause Time 5 s	No. of Repetition 4
Viewing Angle 45°	Model Male Adult	
Category Lips	Equipment none	

Procedures

Exercise



Go to “Trainee List” for verifying the assignment

The screenshot shows the Oroeasy app interface. The top navigation bar includes 'Welcome XXX', 'Exercise', and 'Lips'. The left sidebar contains the following menu items: 'Training Collection', 'Trainee List' (circled in red), 'Language', 'About Oroeasy', 'App Policy', and 'Logout'. The main content area shows a table with columns for 'Pause Time' and 'No. of Repetition'. Four callout boxes provide details for the menu items:

- Training Collection**
 - Contain the training template
- Trainee List**
 - Contain the trainees that are under your care
- Language**
 - Switch among 3 language i.e. English/Traditional Chinese/Simplified Chinese
- Logout**
 - return to the main page



Click “training set” to read the training of the trainee

☰ Trainee List Active ▾

 Reference Code CUHK 001

Registration Status
Completed

 Detail  Training Set

 Reference Code AAA

Registration Status
In Progress

 Detail  Training Set

Note: In the second entry, the "Training Set" button is circled in purple and has a blue arrow pointing to it.



Click “training set” to read the training of the trainee

☰ Trainee List Active ▾

 Reference Code CUHK 001

Registration Status
Completed

 Detail  Training Set

 Reference Code AAA

Registration Status
In Progress

 Detail  Training Set

Note: In the second entry, the "Training Set" button is circled in purple and has a purple arrow pointing to it.



The training info is shown in the panel. The training history has been recorded. Click “training record” to list it.

Lips Strengthening

Start Date: 2021-04-22 End Date: 2021-04-30

Duration: 3 min No. of Exercise: 3

Category: Lips

Progress: 1 / 3

Description: 1. To strengthen the muscles of the lips

Buttons: Delete, Edit, Training Record

Lips Strengthening

Start Date
2021-04-22

End Date
2021-04-30

Progress
 1 / 3

Duration
3 min

No. of Exercise
3

Category
Lips

Description
 1. To strengthen the muscles of the lips

Update

Done
2021-04-20 12:24 - ID:1

Training Record



香港中文大學
The Chinese University of Hong Kong



香港中文大學醫學院
Faculty of Medicine
The Chinese University of Hong Kong



Thank you !